

ClearKlens Cleansinald SC

Sterile Alkylamine surface disinfectant

VH9

Description

ClearKlens Cleansinald SC is a concentrated surface active disinfectant for the Pharmaceutical industry.

Key properties

ClearKlens Cleansinald SC is used for the cleaning and disinfection of walls, floors and work surfaces and in the pharmaceutical industry

ClearKlens Cleansinald SC is a 25ml single shot double bagged disinfectant, which has been filtered to 0.2µm and sterile filled

ClearKlens Cleansinald SC has a broad spectrum of activity

ClearKlens Cleansinald SC can be applied as a spray, wiping or wet sweeping without rinse. Small items or dismantled equipment may be immersed in a

ClearKlens Cleansinald SC solution.

Benefits

Broad spectrum of activity

0.2µm filtered

Double bagged for easy transfer into the cleanroom

Delivered with conformity certificates

Use instructions

ClearKlens Cleansinald SC is used at 0.5% v/v aqueous solution (empty the contents of a 25ml bottle into 5 litres of water).

Technical data

Appearance: clear, slightly yellow liquid

Relative density at 20°C: 1.053

pH (1% solution at 20°C): 9.9

Chemical Oxygen Demand (COD): 563 g O₂ / kg

Nitrogen content (N): 14.5 g/kg

Phosphorus content (P): none

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Store in original closed containers away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet.

Product compatibility

Application, when used as directed, will not effect materials normally encountered in the aforementioned industries.



ClearKlens Cleansinald SC

Microbiological data

ClearKlens Cleansinald SC has been extensively tested using CEN 1040 bactericidal (20°C, 5 minutes), CEN 1276 bactericidal in clean conditions (20°C, 5 minutes) and CEN 1650 fungicidal Activity on *Candida albicans* at 20°C in 15 minutes. Full microbiological data is available in the technical file.

Diversey Europe, Middle East, Africa

P.O. Box 40441, 3504 AE Utrecht

Maarssebroeksedijk 2, 3542 DN

Utrecht

The Netherlands

Tel +31 (0)30 247 69 12

Fax +31 (0)30 247 64 93

www.diverse.com