

## Boma Absorb - 330 ml

### GENERAL

Deodorizing, coagulating absorbent for liquids with acid neutralization.

- Visual neutralization indicator for acids
- Absorbs 100 to 120 times its own volume
- Prevents contact with the spilled liquid
- Can be easily removed with a brush or shovel after solidification
- Particularly suitable for cleaning up hazardous substances or substances with unpleasant odors
- Compact packaging with handy sprinkle cap



### USAGE

- Sprinkle Absorb over the liquids or stains you want to absorb.
- Let the product soak in until it swells and wait until it forms a gel (almost immediate action). For acidic products, wait until the color turns green.
- Remove the gel with a shovel or brush.

Note: on light floor coverings, it is recommended to test the product beforehand.

**Boma Absorb - 330 ml****TECHNICAL SPECIFICATIONS**

Afmetingen product (HxØ) 12.5 x 7 cm

**LOGISTIC DATA****Product**

Article	220019
Unit	Jar
Gross weight	291 g
Net weight	265 g
EAN code	5420015515552

**Packing**

Quantity per packing	6
Dimensions	21.5 x 14 x 15 cm
Gross weight	1833 g

**Pallet**

Quantity per pallet	1848
---------------------	------